

MPOG App Suite Release - December 2024


How to apply the upgrade:

The upgrade package will be sent to each site's IT contact.

- For sites using Desktop Virtualization (e.g. Citrix), your IT team will upgrade the App Suite.
- For sites where the App Suite is installed on PCs, your IT team will distribute the installer
 to all individuals after the database upgrade has been applied. Once you have received
 the installer from your IT team, follow the install instructions <u>here</u>.

DATA TYPES

 Patient Attributes, allowing for greater detail in patient demographic information, is now available. This is optional and may require a supplemental extract. You will not see data in these categories unless your site opts to extract this new data type for submission to MPOG.



DATA DIAGNOSTICS

Added new diagnostic for Anesthesia Assistants



- Retired old individual airway diagnostics in favor of a single combined airway diagnostic
- Updated various diagnostics with new lab unit concepts



CASE VIEWER

- Updated concepts in H&P section
- Fixed bug where "Case Linked" column in the Pro Fee Diagnosis section under the Administrative tab did not display correctly
- Propofol infusion in the case chart are now converted to display in mg instead of mcg
 PROPOFOL
 50 mg
- Updated concept presets in the case chart and global search sections
- Added Microbiology tab. This is optional and may require a supplemental extract. You
 will not see microbiology data unless your site opts to extract this new data type for
 submission to MPOG.



Months in listed dates are now displayed using words (e.g. January 01, 2025)



CASE VALIDATION

Added new concepts to glucose question

NSQIP IMPORT

Removing header metadata has been increased to 50 rows


STS IMPORT

Fixed bugs related to changes in the STS harvest file




VARIABLE MAPPING

- Examine button now includes a "fast mode" setting for sites with large variable counts.
 - Results are random instead of being most recent.



- Updated premappings
- Lab units of measure are now available for mapping



LOCATION MAPPING

Added new "Sub-Unit" location type



TRANSFER TO MPOG CENTRAL

- Upload protocol has been revised. Out of memory issues should be less of an issue.
- Users no longer need to select an upload destination for uploading to TEST vs. PROD

