


## Content Sync


- Suite now checks on startup if any new concepts are available. Notification ball appears if any are available



- Sync now acquires a lock to prevent multiple synchronizations happening at the same time

## Diagnostics

- For some diagnostics, the months of March/April/May 2020 are now excluded from diagnostic pass/fail consideration due to COVID-19's impact on case volume



### Percentage of Outpatient Cases with a LMA Note

Open case list for

- Added Microbiology diagnostic

#### Percentage of Cases with Microbiology Labs

Priority: [Low Priority](#)

Diagnostic Executed On: 2020-09-14 10:05

Open case list for selected month  
(Select a point to enable case listing)

##### Description

Percentage of cases with any microbiology labs. Check extract or mapping if inconsistent with your institution's practices

##### Concepts Used:

###### Microbiology Labs

|      |                                   |
|------|-----------------------------------|
| 3406 | Micro - Culture - Aerobic         |
| 3409 | Micro - Gram Stain                |
| 3418 | Micro - Acid Fast Stain (AFB)     |
| 5155 | Micro - Culture - Aerobic (blood) |
| 5156 | Micro - Culture - AFB (blood)     |