

MPOG App Suite Release – November 2025

How to apply the upgrade:

The upgrade package will be sent to each site's IT contact.

- For sites using Desktop Virtualization (e.g. Citrix), your IT team will upgrade the App Suite.
- For sites where the App Suite is installed on PCs, your IT team will distribute the installer to all individuals after the database upgrade has been applied.
 - Once you have received the installer from your IT team, follow the installation instructions here.
- After applying the upgrade, please confirm that you are using the current version by clicking 'About' in the upper right corner. If your version does not match the below screenshot, please contact support@mpog.zendesk.com.



INFRASTRUCTURE

- Updated labels for CPT codes submitted as part of the hospital discharge diagnoses file (and currently labeled as a discharge billing code) to be listed as "unspecified anesthesia professional fee"
- Intraop, preop, and outcomes notes text values now no longer are truncated to 255 characters
- Patient table now includes dedicated columns for race and ethnicity
- Fixed a bug where registry data was not included in the upload

TRANSFER TO MPOG CENTRAL

- Variable mappings are now included in uploading to MPOG Central
- BRI (Blinded Record Index) has been renamed to PPRL (Privacy Protecting Record Linking) in the interface




CASE VIEWER

 Record search box has been made larger. Also includes text wrapping and 'X' button to clear the search



Floor/ICU documentation is now available in the Table View tab



- Logic for cardiopulmonary bypass timing window has been updated to more closely match the MPOG phenotype
- Record search presets now include the 'mpog:' prefix to help narrow search results to MPOG concepts.



Fixed a bug where the Microbiology tab would display improperly and would have part
of the previous tab still showing.

CASE VALIDATION


Infusion rates now round to two decimal places





VARIABLE MAPPING

- Fixed a crash bug caused by multiple instance IDs using the same organization ID
- Added variable details for the Patient Attributes data type under the "Examine" button



DATA DIAGNOSTICS


- Fixed a bug where certain lab diagnostics would not refresh
- Fixed a bug where the diagnostic graph would not properly update after selecting a Location
- The Fluid Unit Mapping diagnostic previously assessed unit mapping for both inputs and outputs. This could cause the graph results to show 'Poor' performance due to some outputs not having an associated unit. This diagnostic has been updated to no longer include outputs.

IMPORT MANAGER ASSISTANT

Log tabs now default to past one month to improve performance



Overview tab now automatically scrolls to the most recent month





CONTENT SYNCHRONIZATION

• New MPOG concept types are now supported for syncing. This fixes the issue where sites could not successfully run Content Synchronization after an upgrade release until the most recent upgrade was applied at the site.