

## **TIVA Tips**

# Please Scan this QR code and tell us what you think!



### **Avoid Distressing Awareness**

- Check your equipment (e.g. IV is running, infusion pumps are working and correctly programmed).
- Provide adequate analgesia with an opioid infusion (e.g. remifentanil 0.05-0.1 mcg/kg/min) or non-opioid infusion (e.g. dexmedetomidine 0.2-1 mcg/kg/hr)
- Use EEG monitor to guide pharmacodynamic endpoint (sufficient hypnotic effect).
- Target alpha/theta spindles and delta waves on EEG trace or red train tracks in delta and alpha frequency bands on the spectrogram (depending on which monitor is being used).
- Do not rely exclusively on processed EEG index.
- Use pharmacokinetic modeling (e.g. stanpumpR.io) to guide sufficient hypnotic dosing (effect-site concentration).
- Consider targeting <u>brief</u> deep anesthesia (e.g. periods of EEG suppression) for intense stimuli (e.g. intubation, incision)
- Avoid excessive NM-blockade (e.g. maintain 2 twitches on TOF).
- Reverse NM-blockade prior to discontinuing propofol at the end of the case.
- Ensure sufficient analgesia (e.g. sufentanil, fentanyl, methadone, hydromorphone, ketorolac, ketamine) is on board at emergence.

### **Avoid Excessive Hypnosis**

- Use EEG monitor to guide pharmacodynamic endpoint (not excessive hypnotic effect, e.g. EEG suppression).
- Target alpha/theta spindles and delta waves on EEG trace or red train tracks in delta and alpha frequency bands on the spectrogram (depending on which monitor is being used).
- Do not rely exclusively on processed EEG index.
- Use pharmacokinetic modeling (e.g. stanpumpR.io) to avoid excessive hypnotic dosing (effect-site concentration).
- Decrease propofol infusion rate throughout the case (guided by PK modeling) to avoid excessive accumulation.
- Consider decreasing propofol concentration rate in the face of EEG suppression.

#### **Avoid Prolonged Emergence**

- Use EEG monitor to guide down-titration of propofol towards the end of the case, while maintaining sufficient hypnotic effect.
- Do not rely exclusively on processed EEG index.
- Use pharmacokinetic modeling (e.g. stanpumpR.io) to guide safe down-titration of propofol towards the end of the case (target lower effect site concentration) to avoid excessive accumulation.
- Provide adequate analgesia (e.g. remifentanil, sufentanil, fentanyl, methadone, hydromorphone) to allow minimization of propofol.
- Discontinue propofol early while continuing analgesic administration (e.g. remifentanil 0.2 mcg/kg/min) towards the end of the case after reversal of NM-blockade.